



# DualGT

## The Deluxe and Powerful Mattress Replacement System

With the Anti-bottoming Out Sensor providing exceptional protection and interchangeable air cells for total and constant pressure redistribution, DualGT devoted extra attention to the sacrum area, one of the most common sites that develop pressure ulcers. DualGT offers multiple comfort level configurations that deliver optimal internal pressure according to each patient's sensation to create maximum comfort.



### Pressure Indication

The visual indicator demonstrates the optimal pressure level for the patient



### Intelligent Pressure Sensing Technology (IPS)

Responds to patient movement on the mattress by automatically adjusting internal pressure



### Comfort Setting

Comfort adjustment options provide the optimal pressure according to each patient's sensation; there's no need to input weight



### Dual Therapy Modes

Provides Dynamic and Static therapy functions to meet individual patients' needs



### Multiple Cycle Time

Choose 5, 10, or 15 minute cycle times for alternating pressure therapy



### Max Inflate

This feature rapidly brings the mattress to maximum steady pressure, allowing the caregiver to perform the daily nursing procedures

## The DualGT Mattress

The modular air cells of DualGT are independently divided into different zones yet fully interconnected. The Head Zone is always at constant low pressure, providing a pillow-like support function, while in the Torso and Leg Zones, alternating pressure supports the chest and leg areas.

### The Anti-Bottoming Out Sensor

- Extra protection for the sacrum area

### Easily Exchangeable Air Cells

- Detachable air cells to exchange when necessary

### Bridging Cell

- The special replaceable 5" cell can provide constant pressure relief for localized healing.

### Square Air Cell

- Without gaps between air cells, the mattress can deliver maximum support and pressure redistribution and additional safety for getting in and out of the mattress

### Excellent Hygiene Control

- The mattress is completely covered with a durable and washable coverlet, with an easy wipe-clean, hospital-grade mattress base

### Instant CPR

- Rapid delating the entire mattress within seconds in case of emergency

### Effortless Portability Design

- The mattress is thoughtfully designed for easy carrying, featuring color-coded straps for convenient identification when it is folded.



Technical Data	Power Unit Model	REF : SR361
	Power Input	220-240Vac 50 / 60Hz 0.2A
	Power Consumption	Max : 15W
	Weight	<ul style="list-style-type: none"> <li>Power Unit : 2.8Kg ( 6.2lb )</li> <li>Mattress Set : 10kg, 12.5kg ( 22lb, 28lb )</li> </ul>
	Dimension	Power Unit ( W x H x D ) : 12" x 9" x 5" ( 29.5 x 22.5 x 12cm ) Mattress Set ( W x L x H ) : <ul style="list-style-type: none"> <li>32" x 80" x 8" ( 82 x 200 x 21cm )</li> <li>36" x 80" x 8" ( 90 x 200 x 21cm )</li> <li>Other size upon request</li> </ul>
	Agency Approval	<ul style="list-style-type: none"> <li>IEC/EN 60601-1</li> <li>ANSI AAMI ES60601-1</li> <li>AS/NZS IEC60601.1</li> </ul>
	Operating Condition	<ul style="list-style-type: none"> <li>Temperature : 5°C to 40°C ( 41°F to 104°F )</li> <li>Humidity : 15% to 90%</li> </ul>
	Transportation & Storage Condition	<ul style="list-style-type: none"> <li>Temperature : -25°C to 70°C ( -13°F to 158°F )</li> <li>Humidity : 0% to 90%</li> </ul>
Weight Capacity	250kg ( 550lb )	

This Power Unit is not AP/APG protected

[www.carilexmedical.com](http://www.carilexmedical.com)